

Safe Stride Immersive Therapy

Empowering Recovery with AI-Driven Immersive Therapy

Case Study: Virtual Reality for Balance Rehabilitation in Patients with Vestibular Disorders:

Introduction:

Vestibular disorders affect the inner ear and brain, leading to balance and spatial orientation issues. These conditions, such as benign paroxysmal positional vertigo (BPPV), Meniere's disease, and vestibular neuritis, can cause dizziness, vertigo, and imbalance, significantly impacting daily life. Rehabilitation exercises, traditionally administered through physical therapy, have been effective in improving symptoms and decreasing disability. However, Virtual Reality (VR) technology is emerging as a new tool to augment these rehabilitation efforts.

Background on Vestibular Disorders

The vestibular system, located in the inner ear, is responsible for detecting motion, spatial orientation, and maintaining balance. Damage to this system leads to vestibular disorders, which manifest as:

- **Dizziness** or feeling off-balance
- **Vertigo** (a sensation of spinning)
- **Nausea**
- **Difficulty with coordination** or walking

Vestibular rehabilitation therapy (VRT) is a set of exercises designed to improve gaze stability, balance, and motion sensitivity. These exercises aim to help patients adapt and compensate for vestibular dysfunction. However, some patients may find conventional therapy repetitive, boring, or hard to adhere to because of fear of falling or muscle weakness. This is where VR offers an innovative alternative.

Use of Virtual Reality in Vestibular Rehabilitation

VR technology creates immersive environments that can simulate real-world scenarios, allowing for interactive and adaptable rehabilitation experiences. By simulating various environments and movements, VR can challenge patients' balance in a controlled manner.

Key Advantages of VR in Rehabilitation:

1. **Engagement and Motivation:** The immersive nature of VR can make rehabilitation more engaging. Gamified environments or tasks motivate patients to complete exercises with higher adherence rates.
2. **Controlled Environment:** VR allows therapists to manipulate environmental variables (e.g., movement, speed, or visual stimuli), creating a customizable experience that can be gradually intensified.
3. **Safe Simulation of Real-World Scenarios:** Patients can practice tasks that may be dangerous or uncomfortable in real life, such as walking through busy streets, navigating uneven surfaces, or riding in vehicles.
4. **Real-Time Feedback:** Patients receive immediate feedback on their performance, helping them adjust and improve their movements during exercises.
5. **Remote Access:** VR platforms can be used for remote therapy, allowing patients to continue rehabilitation at home under the guidance of a therapist.

Case Presentation:

Patient Profile:

- **Age:** 55-year-old female
- **Diagnosis:** Chronic vestibular neuritis resulting in imbalance and frequent dizziness
- **Symptoms:** Difficulty walking in crowded spaces, vertigo when turning her head quickly, and occasional falls.
- **Previous Treatment:** Traditional vestibular rehabilitation therapy (VRT) for 6 months with moderate improvement but limited progress in dynamic environments.

Intervention:

- The patient was introduced to an 8-week VR-based rehabilitation program, supplementing her traditional VRT. The VR system involved using a head-mounted display (HMD) and balance board with exercises designed to target gaze stabilization, postural control, and balance in dynamic situations.

Rehabilitation Protocol:

1. **Week 1-2:** The initial phase focused on simple balance exercises in a static virtual environment, where the patient stood on the balance board and focused on visual targets that moved slowly.
2. **Week 3-4:** The level of activity and difficulty increased with patient performing leg endurance and coordination activities in a VR exercise platform while actively moving and following instructions.
3. **Week 5-6:** The complexity of the program increased with adding tasks of daily living focusing on precision, coordination, focus and muscle activity.
4. **Week 7-8:** Progressed to a more dynamic activity such as walking from the City center to the neighborhood park allowing her to experience real life situations allowing her to turn

her head, change her speed and also enjoy the scenery around her, Allowing her to exercise her gaze and improve focus.

5. **End of Week 8:** The final stage introduced head movements and rotational exercises to desensitize the vestibular system to fast head motions. The patient "walked" through a virtual forest, turning her head as she looked for specific objects.

Outcome Measures:

- **Objective Balance Measures:** The patient's balance was evaluated using the Berg Balance Scale (BBS) and timed walking tests. Pre-intervention BBS score was 40/56, indicating moderate impairment.
- **Patient-Reported Outcomes:** The Dizziness Handicap Inventory (DHI) and Activities-specific Balance Confidence (ABC) Scale were used to assess perceived disability and confidence in daily activities.

Results:

- After 8 weeks, the patient's BBS score improved to 51/56, reflecting significantly improved balance. Timed walking tests also showed improved lower extremity strength and endurance in navigating dynamic environments.
- On the DHI, the patient's score decreased by 30%, indicating reduced dizziness and better management of symptoms.
- The patient reported feeling more confident walking in crowded spaces and experienced fewer episodes of vertigo when turning her head.

Discussion:

The use of VR in vestibular rehabilitation proved beneficial in this case. The immersive nature of VR enhanced the patient's motivation and engagement, which may have contributed to better adherence compared to traditional VRT. VR's ability to simulate real-world environments allowed the patient to practice and improve her balance in dynamic settings, something difficult to replicate in a clinic.

While VR is not intended to replace traditional therapy, it offers a valuable complement, especially for patients struggling with specific aspects of balance or those seeking more engaging therapy options. The success of this intervention demonstrates VR's potential for broader use in vestibular rehabilitation.

Limitations and Considerations

While the case presented favorable outcomes, there are some considerations to note:

- **Accessibility:** Not all patients have access to VR technology, and the cost can be prohibitive.

- **Simulation Sickness:** Some patients may experience dizziness or nausea from VR use, potentially exacerbating vestibular symptoms. This is called Cybersickness. The best way to counter that is slow progression into the program.
- **Need for Tailored Programs:** Rehabilitation programs should be carefully tailored to each patient's needs, considering their specific symptoms and limitations.

Conclusion:

This case study highlights the promising role of Virtual Reality in the rehabilitation of patients with vestibular disorders. By creating engaging, customizable, and safe environments for balance training, VR offers a novel tool that can improve the outcomes of traditional vestibular rehabilitation. Further research is warranted to explore its long-term benefits and the scalability of VR programs for broader clinical use.

VR is the biggest hope on the world of neurotherapy, neuro-rehab and prevention of falls in our aging population!

Dr. Smith Ahluwalia

Dr. Vik Ahluwalia